

RACK / BOX LABEL / STRAW LABELS




Diagram	Order Reference	Label Width A (mm)	Label Height B (mm)	Write-On Area C (mm)	D (mm)	Labels/roll	Printer Ribbon	Special Properties
▲	BPTLAB-35-427YL*	20.00	20.00	-	-	1000	R-4300	F V
◐	EPT-656-461	30.00	25.40	9.53	-	1000	R-4900	F N A H S L V
◐	EPT-627-461	27.94	6.35	22.00	-	1000	R-4900	F N A H S L

* these labels are also suitable for Cryotheque sample storage software

** AW= All White / no transparent area