

## WELL PLATE LABELS




Diagram	Order Reference	Label Width A (mm)	Label Height B (mm)	C (mm)	D (mm)	Labels/roll	Printer Ribbon	Special Properties
▲	BPT-1-8423	19.05	6.35	-	-	2500	R-7962	F A SPG
▲	BPT-6-8423	38.10	19.05	-	-	2500	R-7962	F A PG
▲	BPT-635-488	48.26	6.35	-	-	2500	R-4300	F A PG
▲	BPT-629-498	38.10	19.05	-	-	2500	R-4300	F R PG