

# LAB PAL™ Labels

Order References			Material	Application
9.5 mm wide	12.7 mm wide	19.1 mm wide		
PAL-375-461	PAL-500-461	PAL-750-461	Polyester B-461 / White*	<b>F</b> <b>N</b> <b>A</b> <b>H</b> <b>L</b> <b>V</b> <b>P</b> <b>G</b>
PAL-375-499	PAL-500-499	PAL-750-499	Nylon Cloth B-499 / White	<b>F</b> <b>N</b> <b>A</b> <b>V</b> <b>P</b> <b>G</b>
PAL-375-488	PAL-500-488	PAL-750-488	Polyester B-488 / White	<b>F</b> <b>A</b> <b>S</b> <b>P</b> <b>G</b>
PAL-375-423	PAL-500-423	PAL-750-423	Polyester B-423 / White	<b>F</b> <b>A</b> <b>P</b> <b>G</b>
PAL-375-430	PAL-500-430	PAL-750-430	Polyester B-430 / Clear	<b>A</b> <b>P</b> <b>G</b>
PAL-375-439	PAL-500-439	PAL-750-439	Vinyl B-439 / White**	<b>F</b> <b>P</b> <b>G</b>

- F** Freezer (-70°C)
- N** Liquid Nitrogen (-196°C)
- A** Autoclave
- H** Hot Water Baths (100°C)
- L** Low-Profile (extremely thin)
- V** Vial Identification
- S** Slide Identification
- P** Plate Identification
- G** General Identification

*B-499 cartridges are 5.8 m long, all others are 6.0 m long*

\* B-461 is acceptable for use in liquid nitrogen only when wrapped upon itself.

\*\* B-439 vinyl is available in Blue, Green, Orange, Red and Yellow. When ordering, please add BL, GR, OR, RD, YL after part number.