



CENTRIFUGE / EPPENDORF / PCR TUBE LABELS

Diagram	Order Reference	Label Width A (mm)	Label Height B (mm)	Write-On Area/ Diameter C (mm)	Labels/ Cartridge	Special Properties
0.2 ml PCR TUBE						
	XSL-108-461	4.75	38.10	17.78	500	FNAHSL VS
0.5 TO 0.7 ml TUBE						
	XSL-73-461	9.53	45.72	19.05	225	FNAHSL VS
	XSL-72-461	12.70	44.45	19.05	225	FNAHSL VS
	XSL-123-461	9.53	28.58	19.05	500	FNAHSL VS
	XSL-122-461	12.70	28.58	19.05	500	FNAHSL VS
	XC-375-461	9.1 m	9.53	-	9.1 m roll	FNAHSL VS
	X-112-490	31.75	9.53	-	500	FNAH L V PG
	X-130-499	20.96	9.53	-	275	FNAH L V PG
	X-118-499	25.40	9.53	9.53	300	FNAH L V PG
	X-120-499	25.40	12.70	9.53	300	FNAH L V PG
1.0 TO 8.0 ml TUBE						
	XSL-74-461	12.70	55.88	19.05	175	FNAHSL VS
	XSL-86-461	15.24	66.68	25.40	150	FNAHSL VS
	XSL-124-461	12.70	41.91	19.05	225	FNAHSL VS
	XSL-126-461	15.24	45.72	25.40	225	FNAHSL VS
	XC-500-461	9.1 m	12.70	-	9.1 m roll	FNAHSL VS
	X-124-490	41.91	12.70	-	225	FNAH L V PG
	X-131-499	25.40	12.70	-	225	FNAH L V PG
	X-118-499	25.40	9.53	9.53	300	FNAH L V PG
	X-120-499	25.40	12.70	9.53	300	FNAH L V PG
	X-117-499	25.40	12.70	11.18	300	FNAH L V PG