


- F** Freezer environment
- N** Liquid nitrogen environment
- A** Autoclave environment
- H** Hot water (100°C) bath
- S** Self-laminating label
- L** Low-profile (extremely thin) label
- R** Removable

- V** Vial identification
- S** Slide identification
- P** Plate identification
- G** General identification



CONICAL / LARGE TUBE LABELS

| Diagram | Order Reference | Label Width A (mm) | Label Height B (mm) | Write-On Area/ C (mm) | Labels/ Cartridge | Special Properties |
|--------------------------------|-----------------|-----------------------|------------------------|--------------------------|----------------------|-----------------------|
| 14.0 – 15.0 ml CONICALS | | | | | | |
| | XSL-119-461 | 38.10 | 95.25 | 31.75 | 100 | FNAHSL VS |
| | XSL-129-461 | 38.10 | 57.15 | 31.75 | 100 | FNAHSL VS |
| | XC-1000-461 | 9.1 m | 25.40 | – | 9.1 m roll | FNAHSL VS |
| | XC-1250-461 | 9.1 m | 31.75 | – | 9.1 m roll | FNAHSL VS |
| | X-113-490 | 38.10 | 55.88 | – | 175 | FNAH L V PG |
| | X-132-499 | 31.75 | 38.10 | – | 200 | FNA V PG |
| 50.0 ml CONICALS | | | | | | |
| | XSL-119-461 | 38.10 | 95.25 | 31.75 | 100 | FNAHSL VS |
| | XC-1000-461 | 9.1 m | 25.40 | – | 9.1 m roll | FNAHSL VS |
| | XC-1250-461 | 9.1 m | 31.75 | – | 9.1 m roll | FNAHSL VS |
| | X-114-490 | 38.10 | 95.25 | – | 100 | FNAH L V PG |
| | X-78-499 | 25.40 | 48.26 | – | 125 | FNA V PG |