


- F** Freezer environment
- N** Liquid nitrogen environment
- A** Autoclave environment
- H** Hot water (100°C) bath
- S** Self-laminating label
- L** Low-profile (extremely thin) label
- R** Removable

- V** Vial identification
- S** Slide identification
- P** Plate identification
- G** General identification



## STRAW LABELS

Diagram	Order Reference	Straw Description	Label Width A (mm)	Label Height B (mm)	Write-On Area/ C (mm)	Labels/ Cartridge	Special Properties
	XSL-133-427	Small Straws	44.45	25.40	9.53	175	<b>F</b> <b>N</b> <b>S</b> <b>V</b> <b>S</b> <b>P</b> <b>G</b>
	XSL-134-427	Large Straws	44.45	38.10	12.70	175	<b>F</b> <b>N</b> <b>S</b> <b>V</b> <b>S</b> <b>P</b> <b>G</b>
	X-105-427-AW**	Straw Rods	27.94	6.350	-	500	<b>F</b> <b>N</b> <b>S</b> <b>V</b> <b>S</b> <b>P</b> <b>G</b>

\*\*AW = All White / no transparent area