



VIAL TOP LABELS

Diagram	Order Reference	Label Width A (mm)	Label Height B (mm)	Diameter C (mm)	Labels/ Cartridge	Special Properties
0.6 TO 1.5 ml VIAL / TUBE						
✖	X-82-499	9.53	–	–	300	FNA V PQ
✖	X-82-499-YL-BK*	9.53	–	–	300	FNA V PQ
✖	X-82-499-OR-BK*	9.53	–	–	300	FNA V PQ
✖	X-82-499-RD-BK*	9.53	–	–	300	FNA V PQ
✚✚	X-118-499	25.40	9.53	9.53	300	FNA V PQ
✚✚	X-120-499	25.40	12.70	9.53	300	FNA V PQ
1.0 TO 2.0 ml VIAL / TUBE						
✖	X-83-499	12.70	–	–	300	FNA V PQ
✖	X-83-499-YL-BK*	12.70	–	–	300	FNA V PQ
✖	X-83-499-OR-BK*	12.70	–	–	300	FNA V PQ
✖	X-83-499-RD-BK*	12.70	–	–	300	FNA V PQ
✚✚	X-118-499	25.40	9.53	9.53	300	FNA V PQ
✚✚	X-120-499	25.40	12.70	9.53	300	FNA V PQ
✚✚	X-117-499	25.40	12.70	11.18	300	FNA V PQ

* Legend will print Black on Coloured Background