



## CONICAL / BOTTLE / PETRI-DISH / DEEP WELL PLATE LABELS

Diagram	Order Reference	Label Width A (mm)	Label Height B (mm)	Write-On Area/ C (mm)	Labels/ Package	Printer Ribbon	Special Properties
<b>14.0 to 15.0 ml CONICAL TUBE</b>							
▲	PTL-113-490	38.10	55.88	-	100	R-4310	FNAHLVPG
▶	PTL-77-461	88.90	38.10	31.75	100	R-4310	FNAHLV
<b>50.0 ML CONICAL TUBE / BOTTLES / PETRI-DISH / DEEP WELL PLATE</b>							
▲	PTL-7-423	12.70	12.70	-	500	R-6010	F A PG
▲	PTL-14-423	25.40	4.75	-	750	R-6010	F A PG
▲	PTL-15-423	25.40	6.35	-	750	R-6010	F A PG
▲	PTL-17-423	24.40	12.70	-	500	R-6010	F A PG
▲	PTL-17-430	25.40	12.70	-	500	R-6010	F A PG
▲	PTL-17-449	25.40	12.70	-	500	R-4310	F A R PG
▲	PTL-17-498	25.40	12.70	-	500	R-6210	F RV PG
⊙	PTL-42-439-OR	25.40	15 m	-	Continuous	R-4310	F PG
⊙	PTL-42-439-YL	25.40	15 m	-	Continuous	R-4310	F PG
▲	PTL-29-498	38.10	12.70	-	500	R-6210	F RV PG
▲	PTL-29-499	38.10	12.70	-	500	R-4310	FNA V PG
▲	PTL-30-427-AW**	38.10	19.05	-	250	R-4310	FN V PG
▲	PTL-30-499	38.10	19.05	-	250	R-4310	FNA V PG
▲	PTL-31-423	38.10	25.40	-	250	R-6210	F A PG
▲	PTL-31-498	38.10	25.40	-	250	R-6210	F RV PG
▲	PTL-31-499	38.10	25.40	-	250	R-4310	FNA V PG
▲	PTL-114-490	38.10	95.25	-	100	R-4310	FNAHLVPG
▲	PTL-78-428-AW**	48.26	25.40	-	250	R-4310	FN V PG
▲	PTL-78-488	48.26	25.40	-	250	R-4310	F A VPG
▲	PTL-78-499	48.26	25.40	-	250	R-4310	FNA V PG
▲	PTL-20-423	50.80	25.40	-	100	R-6010	F A PG
▲	PTL-22-423	76.20	25.40	-	100	R-6010	F A PG
⊙	PTL-8-430	12.70	15 m	-	Continuous	R-6010	A PG

\*\*AW = All White / no transparent area